

## IC Plus Intensity Monitoring Ionisation Chamber

*Ion chambers can be used in a variety of experiments, including: X-ray protein crystallography, generic diffraction experiments, and absorption spectroscopy.*

*A complete system requires a detector head, high voltage supply, and signal processing module. FMB Oxford supplies a full range of products in various configurations, along with the relevant cabling, to suit all your application requirements.*

*The detector head consists of a gas tight chamber containing two parallel plates, over which a high voltage is applied to produce a potential gradient. When X-rays pass through the chamber some are absorbed producing electron-ion pairs that are subsequently separated by the gradient and collected at the plates. The small current produced in the plates is then processed by the electronics module to give a quantitative representation of the X-ray intensity.*

The IC Plus ionisation chambers are suitable for use in a number of diffraction and crystallography experiments, including protein crystallography. The IC Plus is available in 4 different electrode lengths to allow optimisation of the detector energy range and fit the physical space constraints of complex end stations.

The shortest electrode lengths give a compact intensity monitor that does not absorb a large fraction of the beam, ideal for protein crystallography experiments where the end station is typically very crowded and scatter should be minimized. The longer electrode lengths yield more accurate measurements at higher energies and might be more suited to certain spectroscopy applications.



### Specifications

|                     |                              |
|---------------------|------------------------------|
| Window aperture     | 10(h) x 50(w) mm             |
| Windows             | 70µm conducting polymer      |
| Body material       | Aluminium alloy              |
| Electrode gaps      | 10, 14, 18mm                 |
| Electrodes          | Gold plated with guard rings |
| Electrode lengths   | 10, 50, 150, 300mm           |
| Body lengths        | 39, 85, 185, 335mm           |
| Working pressure    | 0.7 - 1.3 Bar Abs            |
| Operating potential | Upto 1.7kV                   |

## Ordering information

|   |          |
|---|----------|
| IC Plus 10, electrode length 10mm, dimensions 98(w) x 39(l) x 66(h)mm     | YMC-0010 |
| IC Plus 50, electrode length 50mm, dimensions 102(w) x 85(l) x 66(h)mm    | YMC-0001 |
| IC Plus 150, electrode length 150mm, dimensions 102(w) x 185(l) x 66(h)mm | YMC-0002 |
| IC Plus 300, electrode length 300mm, dimensions 102(w) x 335(l) x 66(h)mm | YMC-0003 |
| Adapter plate, IC Plus/Prime to ISO80                                     | YMC-0008 |
| Adapter plate, IC Plus/Prime to IC Plus/Prime                             | YMC-0021 |

Your local agent is:

**Head Office**  
FMB Oxford  
Unit 1 Ferry Mills  
Osney Mead  
Oxford OX2 0ES  
United Kingdom

Tel +44 (0)1865 320300  
Fax +44 (0)1865 320301

[enquiries@fmb-oxford.com](mailto:enquiries@fmb-oxford.com)



Certificate Number FM27989

[www.fmb-oxford.com](http://www.fmb-oxford.com)